



Intra-examiner repeatability and inter-examiner reproducibility of subjective refraction

Authors : Kelly Woog, PhD ; Jérôme Perderiset Msc ; Damien Gatineau, MD, PhD ; Laure Pichereau, Msc

Objectifs :

- Literature review on repeatability and reproducibility of subjective refraction to examine the allowable range of dispersion between two measurements,
- Bibliographical study on the parameters which can be at the origin of refraction errors between two measurements.
- Clinical study comparing a classic subjective refraction performed by an expert to that obtained by a new device allowing the automation of the eye examination process (SiVIEW Solution).

I. Bibliographical study

Factors that may cause refraction errors

- Fluctuations in the optical quality of the retinal image (unstable tear film or floating bodies),
- The measuring step.
- Low VA: lower sensitivity: higher measurement uncertainty,
- The examiner's method of questioning, the patient's answers, and how the answers are interpreted by the examiner,
- Accommodative physiological fluctuations and accommodative spasms,
- The size of the pupil, conditioning the depth of field of the eye,
- Pathologies (diabetes, some drugs...).

Intra-examiner repeatability

	Freemant and Hodd 1955	Zadnik et al., 1992	Rosenfield and Chiu 1995	Smith et al 2006	Goss & Grosvenor 1996
N patients	17	40	12	Nc	Literary Review
Âge	19-38 y	20-43 y	23-60 y	Nc	Literary Review
Ametropia	Emmetrope to hypermetrope	-10D à +2 D Cyl 0 to -1.50D	Mean Sph = -1.79D Mean Cyl = -0.45D	Nc	Literary Review
Conclusions	Mean difference in absolute value : • 0.13D (S) • 0.080 (C) • 4° (Axis)	Cycloplegic-free 95% LoA : • ± 0,63 D Note: Confusion between repeatability and reproducibility	95% LoA : • ± 0,30 D for ES • ± 0,27 D for S • ± 0,16 D for C • ± 17.1° for A (attention: mean of the SDs and not of the variances: values lower than those actually obtained)	95% LoA : 0.6 D	80% LoA : • 0.25D min 95% LoA : • 0.50D min (S, ES, C)

Inter-examiner reproducibility

	MacKenzie, 2008	Pesudovs et al., 2007	Grein et al., 2014	Goss & Grosvenor 1996	Leinonen et al., 2006	Sloane et al., 1954	French and Jennings, 1974	Perrigin et al., 1982	Bullimore 1998
N patients	1	16	20	Literary Review	22 healthy patients (T=99)	21	12	32	86
N examiners	40	4	6	Literary Review	8	3	12	3	2
Âge	29 y	Nc	28-72 y	Literary Review	26-89 y	14-18 y	Students	Students	11-60 y
Ametropia	Small myopia	Nc	SE between - 6.19 and +1.98 D Cyl between 0.08 and 4.32 D	Literary Review	-6.1 D à +3.8D	Myopic and weak or no astigmatism	Nc	Nc	SE between -7.50 and +2.50D Cyl max : 4.00D
Conclusions	95% LoA (sur SD) : • ±0,549 D (SE) • ±0,47 D (C) Limit of reproductibility – (√2 x 1.96 x SD) : ±0,78 D (SE)	95% LoA : • ± 0,484 D (SE) • 0,345 D (Cyl) • 0,611 D (total vector) Note: no calculation details	95% LoA : • ±0,2 D à ±0,65 D for ES and cylinder	80% LoA : • 0.25D min 95% LoA : • 0.50D min (S, ES, C)	95% LoA (on T) : • ± 0.74D (ES) • ± 0.34D (C) 95% LoA (0.3<VA<0.45) : • ± 1.14D (SE) 95% LoA (VA ≥ 0.7) : • ± 0.51D (ES) Note: confusion between repeatability and reproducibility	Difference ≤ 0.25D In 73% (SE) In 79% (S) In 81% (C) Difference ≤ 0.50D In 97% (SE) In 90% (S) In 99% (C)	Difference ≤ 0.25D In 68% (S) In 85% (C) Best repeatability when ref between -2.00 and +0.50	Difference ≤ 0.25D In 78% (ES) In 88% (S) In 92% (C) Difference ≤ 0.50D In 95% (ES) In 98% (S) In 98% (C) Axe : 65% dans les 5° 87% dans les 10°	95% LoA : SE • -0.90 à +0.65 D • bias between the 2 clinicians (-0.12D) 95% LoA : Cyl • -0.37D à +0.39D • No bias


Conclusion : Limit of agreement (LoA = 1.96 x SD) in 95% of cases : **+/-0.50D** for the ES, the sphere and the cylinder. **Limit of reproductibility** = √2*LoA ≈ **0.71D** (0,78D for MacKenzie).

II. Clinical Study


Purpose: Comparing a new automated eye examination device to a classical subjective refraction

METHOD

- 83 patients (LE)
- from 15 to 78 y (mean 34 y)
- 1 traditional subjective refraction
- Examiner 1
- 1 automated refraction by SiVIEW
- Examiner 2



- Sphere between -7.75D and +9.50D
- Cylinder between 0 and -4,00D
- ES between -8.50D and +8.75D




Description of SiVIEW solution

- Connected and automated eye examination solution Plug & Play
- Easy, intuitive, fast and accurate
- Possibility of delegating tasks while leaving the practitioner at the heart of the final decision thanks to a detailed report



RESULTS




	S	C	SE
LoA	0,38	0,29	0,39
Repro limit	0,54	0,42	0,55

NEAR VISION :

- Tested on 44 people including 31 with a NV recommended by SiVIEW, according to case study and DV refraction.
- Age : mean 23 y (min 17 y - max 31 y)
- 95% of the additions have been confirmed by the expert and accepted by the patient.
- For only 1 person, the addition was not confirmed by the expert.

Accommodating behaviour (NV)



CONCLUSION

- ✓ No clinical difference between the two measurement methods
- ✓ LoA between 0.29 and 0.39D for sphere, cylinder and SE
- ✓ Reproducibility limit between 0.42 and 0.55D

⇒ Lower than the results found in the literature
⇒ Reliable, precise and reproducible

References :

- Bullimore, M. A., Fusaro, R. E., & Adams, C. W. The repeatability of automated and clinician refraction. *Optometry and vision science: official publication of the American Academy of Optometry*, 1998;75(8), 617-622
- Freeman H. Hodd FAB. Comparative analysis of retinoscopic and subjective refraction. *Br J Physiol Opt* 1955;12:8-36.
- French CN, Jennings JAM. Errors in subjective refraction – an exploratory study. *Ophthalmic Optician* 1974;14(16):797-8.805-6.
- Goss DA & Grosvenor T. Reliability of refraction: a literature review. *J Am Optom Assoc* 1996;67: 619-630.
- Grein, H. J., Schmidt, O., et Ritsche, A. Reproducibility of subjective refraction measurement. *Ophthalmologie*, 2014;111:1057-1064.
- Leinonen J, Laakkonen E, Laatikainen L. Repeatability (test-retest variability) of refractive error measurement in clinical settings. *Acta Ophthalmol Scand*. 2006;84(4):532–536.
- MacKenzie GE. Reproducibility of spherocylindrical prescriptions. *Ophthalmic Physiol Opt*. 2008;28(2):143–150.
- Perrigin J, Perrigin O, Grosvenor T. A comparison of clinical refractive data obtained by three examiners. *Am J Optom Physiol Opt* 1982;59(6):515-9.
- Pesudovs K, Parker KE, Cheng H, Applegate RA. The precision of wavefront refraction compared to subjective refraction and autorefraction. *Optom Vis Sci*. 2007;84(5):387–392.
- Rosenfeld M, Chiu NN. Repeatability of subjective and objective refraction. *Optom Vis Sci*. 1995;72(8):577–579.
- Sloane AE, Ounphy EB, Emmons WV, Gallagher JR. A comparison of refraction results on the same individuals. *Am J Ophthalmol* 1954;37(5):696-9.
- Smith, George. Refraction and visual acuity measurements: what are their measurement uncertainties? *Clinical and Experimental Optometry*, 2006 ;89(2):66-72.
- Zadnik K, Mutti DO, Adams AJ. The repeatability of measurement of the ocular components. *Invest Ophthalmol Vis Sci*. 1992;33(7):2325–2333.