

OPTOMETRIST 88AD

OPHTHALMIC EQUIPMENTS



OVERALL SIZES

Height	1305 mm
Width	605 mm
Depth	820 mm
Rotation	45°R-45°L

HEAD REST

Horizontal range	135 mm
Vertical range	160 mm

SEAT

Up range	180 mm
Minimum height	500 mm
Width	460 mm
Depth	440 mm

BACK-REST

Backward	67°
Height	680 mm
Width	450 mm

Ophthalmic chair OPTOMETRIST with electrical elevation, tilting armrest and backrest, adjustable headrest, horizontal manual sliding or optional electric horizontal sliding

FOOT-REST

Back rotation	90°
---------------	-----

ELECTRICAL SPECIFICATIONS

Tension	220/240 V 110/130 V 50/60 Hz
Power	660 VA

WEIGHT

Shipping weight	80 kg
-----------------	-------

