



CHAIR GEMINI 88HE

OPHTHALMIC EQUIPMENTS

Ophthalmic chair GEMINI with electrical elevation, fixed armrest and optional footrest (see accessories art. 88AP)



OVERALL SIZES

Height	1250 mm
Width	575 mm
Depth	715 mm
Rotation	360°

SEAT

Up range	180 mm
Minimum height	528 mm
Width	475 mm
Depth	430 mm

FIXED BACK-REST

Height	720 mm
Width	475 mm

FOOT-REST

Back rotation	90°
---------------	-----

ELECTRICAL SPECIFICATIONS

Tension	220/240 V 110/130 V
Power	50/60 Hz 660 VA

WEIGHT

Shipping weight	60 kg
-----------------	-------