

COMFORT 88FA

OPHTHALMIC EQUIPMENTS



OVERALL SIZES

Height	1350 mm
Width	650 mm
Depth	880 mm
Rotation	45°R-45°L

HEAD REST

Vertical range	100 mm
----------------	--------

SEAT

Horizontal range	65 mm
Up movement	180 mm
Minimum height	525 mm
Width	480 mm
Depth	460 mm

BACK-REST

Backward	85°
Height	900 mm
Width	480 mm

Ophthalmic chair COMFORT with electrical elevation, tilting armrest and backrest, adjustable headrest, automatic electric tilting drive, optional electric horizontal sliding.

ARM-REST

Arm rest	90°
----------	-----

LEG-REST

Forward	85°
---------	-----

ELECTRICAL

SPECIFICATIONS

Tension	220/240 V 110/130 V 50/60 Hz
Power	660 VA

WEIGHT

Shipping weight	105 kg
-----------------	--------