

SL04

LED 5-STEP MAG SLIT LAMP



Available in the UK through Grafton Optical

sales@graftonoptical.com

01923 233980

SL04

LED SLIT LAMP WITH 5-STEP MAG

- 5-step magnification
- LED illumination
 - More than 50,000 hours working time
 - Reduces stress on the eye
 - Provides a bright and sharp light
- Continuous illumination control
- Yellow filter available

Microscope type	Galilean system, Paralell type, Haag-streit type
Magnification chnager	Rotating drum change
Magnification	6x, 10x, 16x, 25x, 40x
Eyepiece	12.5x
Field of view	6x:33mm; 10x:22mm; 16x:14mm; 25x:8.5mm; 40x:5.5mm
Pupil range	52-82mm
Diopter	+7D
Slit width	0mm - 14mm continuously adjustable
Slit length	1mm - 14mm continuously adjustable
Slit length	0-180 degree continuously adjustable
Slit aperture	Continuous from 1-14mm, 14, 10, 5, 3, 1, 0.2mm
Slit vertical tilt	5°, 10°, 15°, 20°
Filter	Heat absorbtion, grey, red free, cobalt blue
Fixation lamp	Red LED3.5V
Light source	LED
Input voltage	100-240V 48-62Hz
Input power	70VA
Dimensions	740 x 640 x 450 mm
Net weight	21kg
Gross weight	24kg

Grafton Optical Company Limited
Unit 7 River Park Industrial Estate,
Billet Lane, Berkhamsted,
Hertfordshire, HP4 1HL
01923 233980

