



Improve the well-being  
of your eyes



November 2023  
ver. 2 - 2023



Thermo-Pulse Technology

Activa is a new business model generator



Available in the UK from



## Relieve tension symptoms and eye swelling

Our cutting-edge device seamlessly adapts to both professional and home environments.

Designed for swift and enduring relief, it employs controlled heat to efficiently clear meibomian gland debris, coupled with a softening touch for optimal secretion management.

This versatile solution addresses swollen and irritated eyelids caused by blepharitis, meibomian gland dysfunctions, styes, and chalazion. Offering a constant 42-degree Celsius heat for 15 minutes, complemented by a micro-vibration massage, our device ensures heightened relief without compromise.

Safety is paramount, with our device engineered to eliminate any risk of burns. In a world where innovation meets precision, this device stands as a beacon of professional-grade ocular care.

Elevate your eye care standards – invest in excellence today.



- Touch screen controller, user-friendly and easy to manage;
- Through a fully automated procedure the device is able to heat and massage lids providing the best comfort and relief;
- The complete procedure **takes only 15 minutes**; it can be offered simply by using a convenient chair without the presence of professionals;
- Relief from the first application, we recommend its use on a weekly basis.

## Activa is a new business model generator

**Increase customer satisfaction** and save time and money. The Activa benefits can also be enjoyed at home.

**Customize the number of sessions available** for the rental period of your customer.

Scan the QR code

Video  
demonstration



Clinical  
studies



## Are you using an IPL in your clinic?

In the realm of Intense Pulsed Light (IPL) therapy, frequently patients do not manifest the full therapeutic effects following the initial applications. Heating and massaging the meibomian glands stands out as a critical adjunct to IPL therapy.

Enter Activa – the strategic inclusion of this advanced element synergistically amplifies the effects of IPL, promising patients an expeditious and tangible elevation in comfort levels. Step into a new era of ocular care where enhanced patient experience aligns seamlessly with cutting-edge interventions.


**Activa supports IPLs treatment efficacy and IPL patient's relief.**




## Wellness for your eyes

At SBM Sistemi we take the well-being of your eyes very seriously.

- This innovative Eye-mask provides a complete heating and automatic massage technology in a single device, which is **fast** and **painless**;
- The service provided will offer better comfort in daily life for lens wearers, stressed eyes caused by heavy use of digital screens and each application will bring additional relief to the patient;
- **Easily sanitized** and no gel or other disposable required.



**Thermo-Pulse  
Technology**



## Activa benefits



PAINLESS AND RELAXING



CUSTOMIZE NUMBER OF SESSIONS  
AVAILABLE FOR THE RENTAL  
PERIOD OF YOUR CUSTOMER



THE TESTED DISINFECTION  
PROCEDURES ENSURES  
PROTECTION FROM INFECTIONS



NO CONSUMABLES NEEDED



IMMEDIATE FEELING OF RELIEF  
AND WELL-BEING FOR YOUR  
CUSTOMER



HEATING AND MASSAGE  
TECHNOLOGY FOR A FAST  
AND PRACTICAL SOLUTION



ACTIVA INCREASES THE  
EFFECTIVENESS OF IPL TREATMENT  
BY GIVING IMMEDIATE RELIEF



NO SPECIALIST REQUIRED TO  
OPERATE THE DEVICE



CUSTOMIZE THE TREATMENT  
TIME AND VIBRATION  
INTENSITY OF THE EQUIPMENT

Available in the UK from Grafton Optical  
sales@graftonoptical.com  
+44 (0)1923 233980