

SL25 SLIT LAMP



CHARACTERISTICS

The optical lens is made of glass which imported from Germany. All the lens are multi-coated, mildewproof and anti-reflection. All the adaptors are universal, and the eye pieces are of the converging type.

SPECIFICATIONS

Type.....parallel
 Magnification..... Five-step magnifications
 Eye Piece.....12.5X
 Total Magnification.....6X(Ø33mm),10X(Ø 22.5mm),16X(Ø14mm),
 (Diameter of Visual Field)25X(Ø 8.8mm),40X(Ø 5.5mm)
 Diopter Adjustment.....-5D~+3D
 Slit Height.....1mm~8mm continuously adjustable
 Slit Width.....0mm~9mm continuously adjustable
 Light Spot Diameter..... Ø 9, Ø 8, Ø 5, Ø 3, Ø 2, Ø 1, Ø 0.2(mm)
 Slit Angle.....0~180°
 Slit Inclination.....5°,10°,15°,20°
 Filter.....heat absorption, grey, redfree, cobalt blue
 Input Voltage.....110V/220V,60/50Hz
 Input Power.....58W
 Illumination Bulb.....12V 30W halogen bulb
 Fixation Bulb.....(red) LED
 Whole Packing Volume.....0.18m³
 G/Weight.....25kgs



CE FDA

EY107 APPLANATION TONOMETER

SPECIFICATIONS

Measurement Range.....0~10.64Kpa
 Light Ring Displacement.....1.53X2=3.06mm
 Diameter of the testing head.....7mm
 Average Error.....±70pa
 Type.....R type / T type
 Three Prism Head as standard accessory



CE FDA




Prism Head

OPTIONAL ACCESSORIES:




1/3" CCD Adapter



Beam Splitter



Teaching Tube



Yellow Filter



Measuring Eyepiece



Photographic Power Box

GRAFTON OPTICAL COMPANY LIMITED

Crown Hall The Crescent Watford WD18 0QW UK

Tel: 01923 233980 Fax: 01923 246570 E-MAIL: sales@graftonoptical.com

www.graftonoptical.com