


# APPLANATION TONOMETER SL-TM



- Product name : APPLANATION TONOMETER SL-TM
- Intended Use : The Applanation Tonometer is used for measurement of inter-ocular pressure.

■ Specification

Tonometer method	Applanation tonometry (pressure is added on a certain area on the cornea surface)
Installation: B-44 B-45 C-83 C-84	SLIT LAMP SL-102/203/500, SLE070 SLIT LAMP SL-203/500, SLE070 SLIT LAMP SL-50 SLIT LAMP SL-45/55
Tonometric range	0~80 mmHg (one scale 2 mmHg)
Tolerance	± 1.5% of the nominal value or ± 0.5 mmHg, whichever is greater
Weight	400g(B-44), 650g (B-45), 660g(C-83/C-84)

**EU Representative**

**FAUSTO GUZZETTI**  
(SHIN-NIPPON TECHNICAL CENTER)

VIA SOLFERINO,5 21057 OLGiate  
OLONA(VA), Italy

**Supplied by**

**AJINOMOTO®**  
AJINOMOTO TRADING, INC.  
SHIN-NIPPON Medical & Ophthalmic Instruments Dept.

East Wing 7F, TFT Building,  
3-1-22 Ariake, Koto-Ku,  
Tokyo 135-8071, JAPAN  
TEL:813-3528-4416  
FAX: 813-3528-4426  
<http://www.shin-nippon.jp>  
<http://www.ajitrade.com>

**Manufactured by**

 Ohira Co., Ltd

629-961 Niiborishinden,  
Minamiuonuma-Shi, Niigata-Ken  
949-7135, JAPAN