

LIGHT WEIGHT FRAME, ROLY K-0350



MADE IN JAPAN

Features

There are 10 types of PD sizes with 52-70mm(2mm each)and 3 types of colors by PD.


PD66 · 68 · 70 GREEN



PD58 · 60 · 62 · 64 BLUE

PD52 · 54 · 56 PINK OR BLUE

The astigmatism axis can be measured by just turning a knob.



 **Grafton Optical**

Available in the UK through Grafton Optical
sales@graftonoptical.com
01923 233980